








INSTRUMENT INDEX - EMS 7A




| | | | | | | | | | | | | | | |
|--------|-------------------------|-----------|-----------|--|--|--|--|--|--|--|---|--|---|--|
| | | | | | | | | | | |   | |   | |
| 0 | ISSUED FOR CONSTRUCTION | | | | | | | | | | DOC NR: 569-DB7A-AIC-310-002 | CLIENT NR: PRD-AIC-LIS-011 | | |
| B | 90% DD ISSUE | | | | | | | | | | | | | |
| A | 30% DD ISSUE | | | | | | | | | | | | | |
| REV | DESCRIPTION | | | | | | | | | | CLIENT: TAKEDA / BAXALTA PROJECT: BURITI EPCMV PROJECT <div>REV.: 0</div> | | | |
| REVIEW | REV. A | REV. B | REV. C | | | | | | | | | | | |
| DATE | 08FEB2021 | 04MAI2021 | 30JUL2021 | | | | | | | | | | | |
| EXEC | MAB | MAB | MAB | | | | | | | | | | | |
| CHECK | MAF | MAF | MAF | | | | | | | | | | | |
| APPROV | MAJ | RSP | RSP | | | | | | | | | | | |

| <div><div><div><div></div><div>TESSLER</div><div>engenharia</div></div><div><div>ips</div></div></div></div> | | | <div><div><div><div></div><div>Takeda</div></div><div><div></div><div>Hemobrás</div><div>Empresa Brasileira de Hemoderivação e Biotecnologia</div></div></div></div> | | TITLE: INSTRUMENT INDEX - EMS 7A | | | | DOC NR.: 569-DB7A-AIC-310-002 | | SHEET.: 3 / 6 |
|--|-----|---------------|--|---------------|---|----------------------|----------------------------|----------------|----------------------------------|--|------------------|
| | | | | | | | | | CLIENT NR.: PRD-AIC-LIS-011 | | REV.: 0 |
| ID | REV | TAG Nr. | SERVICE | P&ID | LOC. (*) | SPECIFICATION SHEETS | ELEMENTARY DRAWING | ROUTES DRAWING | NOTES | | |
| | | | | LINE OR EQUIP | | | WIRING DIAGRAM | | | | |
| 1 | 0 | TIT-A1002-01 | TEMPERATURE TRANSMITTER - PRINTED MATERIALS ROOM A1002 | 7A-I-1-3-11 | F | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-1-3-11 | 1 | | |
| 2 | 0 | MIT-A1002-01 | HUMIDITY TRANSMITTER - PRINTED MATERIALS ROOM A1002 | 7A-I-1-3-11 | C | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-1-3-11 | 1 | | |
| 3 | 0 | TIT-A1004-01 | TEMPERATURE TRANSMITTER - INBOUND STAGING ROOM A1004 | 7A-I-1-3-11 | F | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-1-3-15 | 1 | | |
| 4 | 0 | MIT-A1004-01 | HUMIDITY TRANSMITTER - INBOUND STAGING ROOM A1004 | 7A-I-1-3-11 | C | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-1-3-15 | 1 | | |
| 5 | 0 | TIT-A1006-01 | TEMPERATURE TRANSMITTER - STAGING INSPECTED - ROOM A1006 | 7A-I-1-3-11 | F | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-1-3-15 | 1 | | |
| 6 | 0 | MIT-A1006-01 | HUMIDITY TRANSMITTER - STAGING INSPECTED - ROOM A1006 | 7A-I-1-3-11 | C | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-1-3-15 | 1 | | |
| 7 | 0 | TIT-A1007-01 | TEMPERATURE TRANSMITTER - CONDITIONING STAGING - ROOM A1007 | 7A-I-1-3-11 | F | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-1-3-15 | 1 | | |
| 8 | 0 | MIT-A1007-01 | HUMIDITY TRANSMITTER - CONDITIONING STAGING - ROOM A1007 | 7A-I-1-3-11 | C | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-1-3-15 | 1 | | |
| 9 | 0 | TIT-A1019-01 | TEMPERATURE TRANSMITTER - VISUAL INSPECTION ROOM A1019 | 7A-I-1-3-11 | F | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-1-3-15 | 1 | | |
| 10 | 0 | MIT-A1019-01 | HUMIDITY TRANSMITTER - VISUAL INSPECTION ROOM A1019 | 7A-I-1-3-11 | C | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-1-3-15 | 1 | | |
| 11 | 0 | TIT-A1029-01 | TEMPERATURE TRANSMITTER - PACKAGING ROOM A1029 | 7A-I-1-3-11 | F | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-1-3-15 | 1 | | |
| 12 | 0 | MIT-A1029-01 | HUMIDITY TRANSMITTER - PACKAGING ROOM A1029 | 7A-I-1-3-11 | C | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-1-3-15 | 1 | | |
| 13 | 0 | TIT-A1035-01 | TEMPERATURE TRANSMITTER - 2-8°C STORAGE | 7A-I-1-3-11 | F | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-1-3-15 | 1 | | |
| 14 | 0 | TIT-A1037-01 | TEMPERATURE TRANSMITTER - STAGING ROOM A1037 | 7A-I-1-3-11 | C | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-1-3-15 | 1 | | |
| 15 | 0 | MIT-A1037-01 | HUMIDITY TRANSMITTER - STAGING ROOM A1037 | 7A-I-1-3-11 | F | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-1-3-15 | 1 | | |
| 16 | 0 | TIT-A1040-01 | TEMPERATURE TRANSMITTER - AMBIENT STORAGE ROOM A1040 | 7A-I-1-3-11 | C | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-1-3-15 | 1 | | |
| 17 | 0 | MIT-A1040-01 | HUMIDITY TRANSMITTER - AMBIENT STORAGE ROOM A1040 | 7A-I-1-3-11 | F | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-1-3-15 | 1 | | |
| 18 | 0 | PDIT-A2001-01 | PRESSURE DIFERENTIAL TRANSMITTER - CAPPING/TRAY OFF - ROOM A2001 | 7A-I-1-3-11 | C | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 | | |
| 19 | 0 | TIT-A2001-01 | TEMPERATURE TRANSMITTER - CAPPING/TRAY OFF - ROOM A2001 | 7A-I-1-3-11 | F | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 | | |
| 20 | 0 | MIT-A2001-01 | HUMIDITY TRANSMITTER - CAPPING/TRAY OFF - ROOM A2001 | 7A-I-1-3-11 | C | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 | | |
| 21 | 0 | PDIT-A2002-01 | PRESSURE DIFERENTIAL TRANSMITTER - ASSEPTIC STAGING - ROOM A2002 | 7A-I-1-3-11 | F | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 | | |
| 22 | 0 | TIT-A2002-01 | TEMPERATURE TRANSMITTER - ASSEPTIC STAGING - ROOM A2002 | 7A-I-1-3-11 | C | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 | | |
| 23 | 0 | MIT-A2002-01 | HUMIDITY TRANSMITTER - ASSEPTIC STAGING ROOM A2002 | 7A-I-1-3-11 | F | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 | | |




(*) LOCATION C=FIELD F=IN FRONT OF PAINEL P=OPERATOR CONSOLE I=PLC

|  | | |   | TITLE: INSTRUMENT INDEX - EMS 7A | | | DOC NR.: 569-DB7A-AIC-310-002 | | SHEET.: 4 / 6 |
|---|-----|---------------|--|--|-------------|----------------------|---|----------------|-------------------|
| | | | | | | | CLIENT NR.: PRD-AIC-LIS-011 | | REV.: 0 |
| ID | REV | TAG Nr. | SERVICE | P&ID | LOC. (*) | SPECIFICATION SHEETS | ELEMENTARY DRAWING | ROUTES DRAWING | NOTES |
| | | | | LINE OR EQUIP | | | WIRING DIAGRAM | | |
| 24 | 0 | PDIT-A2023-01 | PRESSURE DIFERENTIAL TRANSMITTER - QC STORAGE - ROOM A2023 | 7A-I-1-3-11 | C | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 |
| 25 | 0 | TIT-A2023-01 | TEMPERATURE TRANSMITTER - QC STORAGE ROOM A2023 | 7A-I-1-3-11 | F | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 |
| 26 | 0 | MIT-A2023-01 | HUMIDITY TRANSMITTER - QC STORAGE ROOM A2023 | 7A-I-1-3-11 | C | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 |
| 27 | 0 | PDIT-A2025-01 | PRESSURE DIFERENTIAL TRANSMITTER- FILLING ROOM A2025 | 7A-I-1-3-11 | F | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 |
| 28 | 0 | TIT-A2025-01 | TEMPERATURE TRANSMITTER- FILLING ROOM A2025 | 7A-I-1-3-11 | C | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 |
| 29 | 0 | MIT-A2025-01 | HUMIDITY TRANSMITTER- FILLING ROOM A2025 | 7A-I-1-3-11 | F | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 |
| 30 | 0 | PDIT-A2036-01 | PRESSURE DIFERENTIAL TRANSMITTER - FORMULATION ROOM - ROOM A2036 | 7A-I-1-3-11 | C | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 |
| 31 | 0 | TIT-A2036-01 | TEMPERATURE TRANSMITTER - FORMULATION ROOM - ROOM 2036 | 7A-I-1-3-11 | F | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 |
| 32 | 0 | MIT-A2036-01 | HUMIDITY TRANSMITTER - FORMULATION ROOM - ROOM A2036 | 7A-I-1-3-11 | C | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 |
| 33 | 0 | PDIT-A2037-01 | PRESSURE DIFERENTIAL TRANSMITTER - BUFFER PREP - ROOM A2037 | 7A-I-1-3-11 | F | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 |
| 34 | 0 | MIT-A2037-01 | HUMIDITY TRANSMITTER - BUFFER PREP ROOM A2037 | 7A-I-1-3-11 | C | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 |
| 35 | 0 | TIT-A2037-01 | TEMPERATURE TRANSMITTER - BUFFER PREP ROOM A2037 | 7A-I-1-3-11 | C | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 |
| 36 | 0 | TIT-A2006-01 | TEMPERATURE TRANSMITTER - MATERIAL STAGING ROOM A2006 | 7A-I-1-3-11 | F | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 |
| 37 | 0 | MIT-A2006-01 | HUMIDITY TRANSMITTER - MATERIAL STAGING ROOM A2006 | 7A-I-1-3-11 | C | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 |
| 38 | 0 | PDIT-A2012-01 | PRESSURE DIFERENTIAL TRANSMITTER- MATERIAL STAGING - ROOM A2012 | 7A-I-1-3-11 | C | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 |
| 39 | 0 | TIT-A2012-01 | TEMPERATURE TRANSMITTER - MATERIAL STAGING ROOM A2012 | 7A-I-1-3-11 | F | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 |
| 40 | 0 | MIT-A2012-01 | HUMIDITY TRANSMITTER - MATERIAL STAGING ROOM A2012 | 7A-I-1-3-11 | C | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 |
| 41 | 0 | TIT-A2014-01 | TEMPERATURE TRANSMITTER - STAGING ROOM A2014 | 7A-I-1-3-11 | C | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 |
| 42 | 0 | MIT-A2014-01 | HUMIDITY TRANSMITTER - STAGING ROOM A2014 | 7A-I-1-3-11 | F | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 |
| 43 | 0 | TIT-A2008-01 | TEMPERATURE TRANSMITTER - STORAGE ROOM A2008 | 7A-I-1-3-11 | C | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 |
| 44 | 0 | PDIT-A2003-01 | PRESSURE DIFERENTIAL TRANSMITTER - CLEAN EQUIPAMENTE - ROOM A2003 | 7A-I-1-3-11 | C | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 |
| 45 | 0 | PDIT-A2004-01 | PRESSURE DIFERENTIAL TRANSMITTER - PRE WASH - ROOM A2004 | 7A-I-1-3-11 | F | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 |
| 46 | 0 | PDIT-A2005-01 | PRESSURE DIFERENTIAL TRANSMITTER - MAL OUT - ROOM A2005 | 7A-I-1-3-11 | C | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 |

(*) LOCATION C=FIELD F=IN FRONT OF PAINEL P=OPERATOR CONSOLE I=PLC

|  | | |   | TITLE: INSTRUMENT INDEX - EMS 7A | | | DOC NR.: 569-DB7A-AIC-310-002 | | SHEET.: 5 / 6 |
|---|-----|---------------|--|--|-------------|----------------------|---|----------------|-------------------|
| | | | | | | | CLIENT NR.: PRD-AIC-LIS-011 | | REV.: 0 |
| ID | REV | TAG Nr. | SERVICE | P&ID | LOC. (*) | SPECIFICATION SHEETS | ELEMENTARY DRAWING | ROUTES DRAWING | NOTES |
| | | | | LINE OR EQUIP | | | WIRING DIAGRAM | | |
| 47 | 0 | PDIT-A2005-02 | PRESSURE DIFERENTIAL TRANSMITTER - MAL OUT - ROOM A2005 | 7A-I-1-3-11 | C | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 |
| 48 | 0 | PDIT-A2009-01 | PRESSURE DIFERENTIAL TRANSMITTER- CIP/SIP - ROOM A2009 | 7A-I-1-3-11 | F | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 |
| 49 | 0 | PDIT-A2017-01 | PRESSURE DIFERENTIAL TRANSMITTER - GOWN DEGOWN - ROOM A2017 | 7A-I-1-3-11 | C | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 |
| 50 | 0 | PDIT-A2017-02 | PRESSURE DIFERENTIAL TRANSMITTER - GOWN DEGOWN - ROOM A2017 | 7A-I-1-3-11 | C | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 |
| 51 | 0 | PDIT-A2026-01 | PRESSURE DIFERENTIAL TRANSMITTER - MAL OUT - ROOM A2026 | 7A-I-1-3-11 | F | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 |
| 52 | 0 | PDIT-A2027-01 | PRESSURE DIFERENTIAL TRANSMITTER - MAL IN - ROOM A2027 | 7A-I-1-3-11 | C | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 |
| 53 | 0 | PDIT-A2027-02 | PRESSURE DIFERENTIAL TRANSMITTER - MAL IN - ROOM A2027 | 7A-I-1-3-11 | C | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 |
| 54 | 0 | PDIT-A2028-01 | PRESSURE DIFERENTIAL TRANSMITTER - JC ROOM A2028 | 7A-I-1-3-11 | F | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 |
| 55 | 0 | PDIT-A2029-01 | PRESSURE DIFERENTIAL TRANSMITTER - CLEAN STORAGE - ROOM A2029 | 7A-I-1-3-11 | C | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 |
| 56 | 0 | PDIT-A2029-02 | PRESSURE DIFERENTIAL TRANSMITTER - CLEAN STORAGE - ROOM A2029 | 7A-I-1-3-11 | C | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 |
| 57 | 0 | PDIT-A2030-01 | PRESSURE DIFERENTIAL TRANSMITTER - AIR LOCK - ROOM A2030 | 7A-I-1-3-11 | F | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 |
| 58 | 0 | PDIT-A2030-02 | PRESSURE DIFERENTIAL TRANSMITTER - AIR LOCK - ROOM A2030 | 7A-I-1-3-11 | C | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 |
| 59 | 0 | PDIT-A2031-01 | PRESSURE DIFERENTIAL TRANSMITTER - PAL OUT - ROOM A2031 | 7A-I-1-3-11 | C | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 |
| 60 | 0 | PDIT-A2031-02 | PRESSURE DIFERENTIAL TRANSMITTER - PAL OUT - ROOM A2031 | 7A-I-1-3-11 | F | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 |
| 61 | 0 | PDIT-A2032-01 | PRESSURE DIFERENTIAL TRANSM - PAL IN ROOM A2032 | 7A-I-1-3-11 | C | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 |
| 62 | 0 | PDIT-A2032-02 | PRESSURE DIFERENTIAL TRANSM - PAL IN ROOM A2032 | 7A-I-1-3-11 | C | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 |
| 63 | 0 | PDIT-A2034-01 | PRESSURE DIFERENTIAL TRANSMITTER - PAL OUT - ROOM A2034 | 7A-I-1-3-11 | F | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 |
| 64 | 0 | PDIT-A2034-02 | PRESSURE DIFERENTIAL TRANSMITTER - PAL OUT - ROOM A2034 | 7A-I-1-3-11 | C | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 |
| 65 | 0 | PDIT-A2035-01 | PRESSURE DIFERENTIAL TRANSMITTER - MAL OUT - ROOM A2035 | 7A-I-1-3-11 | C | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 |
| 66 | 0 | PDIT-A2035-02 | PRESSURE DIFERENTIAL TRANSMITTER - 2 - MAL OUT - ROOM A2035 | 7A-I-1-3-11 | F | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 |
| 67 | 0 | PDIT-A2038-01 | PRESSURE DIFERENTIAL TRANSMITTER - PAL IN - ROOM A2038 | 7A-I-1-3-11 | C | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 |
| 68 | 0 | PDIT-A2038-02 | PRESSURE DIFERENTIAL TRANSMITTER - PAL IN - ROOM A2038 | 7A-I-1-3-11 | C | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 |
| 69 | 0 | PDIT-A2039-01 | PRESSURE DIFERENTIAL TRANSMITTER - MAL IN - ROOM A2039 | 7A-I-1-3-11 | F | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 |

(*) LOCATION C=FIELD F=IN FRONT OF PAINEL P=OPERATOR CONSOLE I=PLC

|  | | |   | TITLE: INSTRUMENT INDEX - EMS 7A | | | DOC NR.: 569-DB7A-AIC-310-002 CLIENT NR.: PRD-AIC-LIS-011 | | SHEET.: 6 / 6 REV.: 0 |
|---|-----|----------------|--|---|-------------|----------------------|--|----------------|--------------------------------|
| ID | REV | TAG Nr. | SERVICE | P&ID | LOC. (*) | SPECIFICATION SHEETS | ELEMENTARY DRAWING | ROUTES DRAWING | NOTES |
| | | | | LINE OR EQUIP | | | WIRING DIAGRAM | | |
| 70 | 0 | PDIT-A2039-02 | PRESSURE DIFERENTIAL TRANSMITTER - MAL IN - ROOM A2039 | 7A-I-1-3-11 | C | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 |
| 71 | 0 | PDIT-A2040-01 | PRESSURE DIFERENTIAL TRANSMITTER - JC ROOM A2040 | 7A-I-1-3-11 | C | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 |
| 72 | 0 | PDIT-A2041-01 | PRESSURE DIFERENTIAL TRANSMITTER - COMPONENT LOAD - ROOM A2041 | 7A-I-1-3-11 | F | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 |
| 73 | 0 | PDIT-A2041-02 | PRESSURE DIFERENTIAL TRANSMITTER - COMPONENT LOAD - ROOM A2041 | 7A-I-1-3-11 | F | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 |
| 74 | 0 | PDIT-A2046-01 | PRESSURE DIFERENTIAL TRANSMITTER - MAL IN - ROOM A2046 | 7A-I-1-3-11 | C | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 |
| 75 | 0 | PDIT-A2025-03 | PRESSURE DIFERENTIAL TRANSMITTER-FILLING ROOM A2025 | 7A-I-1-3-11 | F | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 |
| 76 | 0 | TIT-FS-2001-01 | TEMPERATURE TRANSMITTER FREEZER FS-2001 - ROOM A2049 | 7A-I-1-3-11 | C | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 |
| 77 | 0 | TIT-FS-2002-01 | TEMPERATURE TRANSMITTER FREEZER FS-2002 - ROOM A2014 | 7A-I-1-3-11 | F | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 |
| 78 | 0 | TIT-FR-3801-01 | TEMPERATURE TRANSMITTER REFRIGERATOR FR-3801 - ROOM A2037 | 7A-I-1-3-11 | C | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 |
| 78 | 0 | TIT-FR-3901-01 | TEMPERATURE TRANSMITTER REFRIGERATOR FR-3901 - ROOM A2036 | 7A-I-1-3-11 | C | PRD-AIC-DSH-001 | 7A-I-0-3-11 7A-I-0-3-07 | 7A-I-2-3-25 | 1 |

(*) LOCATION C=FIELD F=IN FRONT OF PAINEL P=OPERATOR CONSOLE I=PLC